i-STAT hs-Tnl CARTRIDGE

Clinical evaluation with i-STAT hs-TnI to aid in the diagnosis of Myocardial Infarction (MI)

This speaker program is sponsored by and on behalf of Abbott



i-STAT hs-TnI Cartridge

INTENDED USE

The i-STAT hs-TnI cartridge with the i-STAT 1 System is intended for use in the in vitro quantification of cardiac troponin I (cTnI) in whole blood or plasma samples in point of care or clinical laboratory settings. The i-STAT hs-TnI cartridge with the i-STAT 1 System is intended to be used as an aid in the diagnosis of myocardial infarction (MI).

Largest POC high-sensitivity Troponin Assay Study to date*

- A multicenter, prospective study was conducted at 28 clinical sites to evaluate the clinical performance of the *i-STAT hs-TnI* cartridge with *the i-STAT 1 System*.*
- Venous whole blood and plasma specimens were collected with lithium heparin anticoagulant from subjects presenting to the emergency department (ED) with signs and symptoms suggestive of acute coronary syndrome (ACS).

*Note: i-STAT clinical study included 3,585 patients, Siemens included 1,089 and Quidel included 422 Sources:

- i-STAT hs-Tnl Instructions for Use (US)
- Siemens Healthineers Atellica VTLi hs-cTnl Instructions for Use (OUS)
- Quidel TriageTrue High Sensitivity Troponin I Instructions for Use (OUS)







i-STAT hs-TnI Cartridge

i-STAT hs-Tnl CARTRIDGE WITH THE i-STAT 1 SYSTEM

Clinical Performance Data

Abbott Point of Care Pivotal Study: CS-2020-0012

CLINICAL PERFORMANCE | WHOLE BLOOD

Sex-Specific Cut-off Analysis (Female: 13 ng/L | Male: 28 ng/L)

SEX	TIME POINT (HOURS)*	N	MI	NON-MI	SENSITIVITY (%)	SPECIFICITY(%)	PPV(%)	NPV(%)
	0 to 1	1870	129	1741	91.47	83.23	28.78	99.25
FFRANIF	>1 to 3	1799	119	1680	96.64	82.14	27.71	99.71
FEMALE	>3 to 6	724	70	654	97.14	77.83	31.92	99.61
	>6	60	16	44	100.00	54.55	44.44	100.00
	0 to 1	1090	130	960	79.23	84.17	40.39	96.77
MALE	>1 to 3	1025	118	907	90.68	83.90	42.29	98.58
IVIALE	>3 to 6	439	69	370	94.20	82.97	50.78	98.71
	>6	47	12	35	91.67	57.14	42.31	95.24

^{*}All timepoints are relative to ED presentation.

CLINICAL PERFORMANCE | WHOLE BLOOD

Overall Cut-off Analysis: (OVERALL: 21 ng/L)

SEX	TIME POINT (HOURS)*	N	MI	NON-MI	SENSITIVITY (%)	SPECIFICITY(%)	PPV(%)	NPV(%)
	0 to 1	1870	129	1741	86.05	89.37	37.50	98.86
CENANI E	>1 to 3	1799	119	1680	92.44	89.70	38.87	99.41
FEMALE	>3 to 6	724	70	654	95.71	85.78	41.88	99.47
	>6	60	16	44	93.75	65.91	50.00	96.67
	0 to 1	1090	130	960	83.08	78.33	34.18	97.16
NANE	>1 to 3	1025	118	907	92.37	77.95	35.28	98.74
MALE	>3 to 6	439	69	370	95.65	74.32	40.99	98.92
	>6	47	12	35	91.67	54.29	40.74	95.00

^{*}All timepoints are relative to ED presentation.

Dr. Cullen is a researcher and practicing emergency medicine physician in Australia. Not all products discussed in this presentation are available for use in the United States.

The views expressed in this presentation belong to the speaker alone and do not necessarily reflect the views of Abbott or other organizations named in the speaker.

The *i-STAT hs-TnI* cartridge has not been evaluated or cleared by the FDA for use in any accelerated diagnostic protocol. The *i-STAT hs-TnI* cartridge is not cleared for use on i-STAT Alinity in the United States, but Dr. Cullen's study was done in Australia, where the *i-STAT hs-TnI cartridge* is cleared for use on the i-STAT Alinity. The *i-STAT hs-TnI cartridge* is only cleared for use on i-STAT 1 in the United States.

For U.S. audience only
For complete product information
visit globalpointofcare.abbott
APOC-25002638.1

HIGH SENSITIVITY POC TROPONIN TESTING



and AMI diagnosis




PROF LOUISE CULLEN

MBBS PhD FACEM FESC

EMERGENCY PHYSICIAN





Disclosures

The views expressed here are solely my own and do not reflect the opinions or positions of any affiliated employer, organization or entity.

Representative roles: National Heart Foundation, International Federation of Clinical Chemistry, ACVC Cardiac Biomarker Committee.

Advisory board/Committee: Abbott Diagnostics, Siemens Healthineers, Quidel/Ortho Diagnostics

Consultant: Abbott Diagnostics, Siemens Healthineers, Quidel/Ortho Diagnostics, Roche

Institutional grant support: Abbott Diagnostics, Siemens Healthineers, Beckman Coulter, Radiometer Pacific, Roche.

IP/Patent: None

Stocks/Bonds: None

Learning Objectives

- Describe the clinical value of using pointof-care high-sensitivity troponin-I (hs-TnI) assays in the early assessment of patients with suspected acute myocardial infarction in rural institutions
- Identify key challenges in transitioning from contemporary troponin to highsensitivity troponin assays
- Discuss how the implementation of hs-Tnl testing at presentation can improve ED workflows and support timely clinical decision-making
- Discuss approaches to integrating hs-Tnl point-of-care testing into existing chest pain pathways and protocols

The views expressed in this presentation belong to the speaker alone and do not necessarily reflect the views of Abbott or other organizations named in the speaker.

OVERVIEW

Emergency medicine and ACS

Rapid risk assessment strategies


Troponin assays
POC HS assays

Implementation of POC HS troponin assays













Suspected ACS

Impact on the whole of ED





Chest pain

Leading cause of ED presentations

•

<15% with AMI




5

AMI DIAGNOSIS

Combination of:

- Clinical assessment
- ECG
- Cardiac troponin (cTn) testing
- Other investigations














ASSESSMENT STRATEGIES



HighSTEACS



0/1 and 0/2h strategies

Patient presents with a suspected NSTEMI and without an indication for immediate invasive angiography Take hs-cTn at 0 h and 1 h/2 h Patients who do not meet High initial hs-cTn the criteria for either of OR OR the other two pathways Increase in 1 h/2 h hs-cTn Low initial hs-cTn and no increase in 1 h/2 h hs-cTn Appropriate management can be determined based on the hs-cTn levels and clinical situation Rule-out pathway Observe pathway Rule-in pathway ESC-

The *i-STAT hs-TnI* cartridge has not been evaluated or cleared by the FDA for use in any accelerated diagnostic protocol.

2023 **ESC** Guidelines for the management of acute coronary syndromes, EHJ (2023) 00, 1-107; https://doi.org/10.1093/eurheartj/e had191

Table \$4 Assay specific cut-off levels in ng/L within the 0 h/1 h and 0 h/2 h algorithms

0 h/1 h algorithm	Very low	Low	No 1 hΔ	High	1 h∆
hs-cTnT (Elecsys; Roche)	<5	<12	<3	≥52	≥5
hs-cTnl (Architect; Abbott)	<4	<5	<2	≥64	≥6
hs-cTnl (Centaur; Siemens)	<3	<6	<3	≥120	≥12
hs-cTnl (Access; Beckman Coulter)	<4	<5	<4	≥50	≥15
hs-cTnl (Clarity; Singulex)	<1	<2	<1	≥30	≥6
hs-cTnl (Vitros; Clinical Diagnostics)	<1	<2	<1	≥40	≥4
hs-cTnl (Pathfast; LSI Medience)	<3	<4	<3	≥90	≥20
hs-cTnl (TriageTrue; Quidel)	<4	<5	<3	≥60	≥8
hs-cTnl (Dimension EXL; Siemens)	<9	<9	<5	≥160	≥100
0 h/2 h algorithm	Very low	Low	No 2 h∆	High	2 h∆
hs-cTnT (Elecsys; Roche)	<5	<14	<4	≥52	≥10
hs-cTnl (Architect; Abbott)	<4	<6	<2	≥64	≥15
hs-cTnl (Centaur; Siemens)	<3	<8	<7	≥120	≥20
hs-cTnl (Access; Beckman Coulter)	<4	<5	<5	≥50	≥20
hs-cTnl (Clarity; Singulex)	<1	TBD	TBD	≥30	TBD
hs-cTnl (Vitros; Clinical Diagnostics)	<1	TBD	TBD	≥40	TBD E
hs-cTnl (Pathfast; LSI Medience)	<3	TBD	TBD	≥90	TBD Ç
hs-cTnl (TriageTrue; Quidel)	<4	TBD	TBD	≥60	TBD ©

Table \$4 Assay specific cut-off levels in ng/L within the 0 h/1 h and 0 h/2 h algorithms

0 h/1 h algorithm	Very low	Low	N o 1 hΔ	High	1 h∆
hs-cTnT (Elecsys; Roche)	<5	<12	<3	≥52	≥5
hs-cTnl (Architect; Abbott)	<4	<5	<2	≥64	≥6
hs-cTnl (Centaur; Siemens)	<3	<6	<3	≥120	≥12
hs-cTnl (Access; Beckman Coulter)	<4	<5	<4	≥50	≥15
hs-cTnl (Clarity; Singulex)	<1	<2	<1	≥30	≥6
hs-cTnl (Vitros; Clinical Diagnostics)	<1	<2	<1	≥40	≥4
hs-cTnl (Pathfast; LSI Medience)	<3	<4	<3	≥90	≥20
hs-cTnl (TriageTrue; Quidel)	<4	<5	<3	≥60	≥8
hs-cTnl (Dimension EXL; Siemens)	<9	<9	<5	≥160	≥100
0 h/2 h algorithm	Very low	Low	No 2 h∆	High	2 h∆
hs-cTnT (Elecsys; Roche)	<5	<14	<4	≥52	≥10
hs-cTnI (Architect; Abbott)	<4	<6	<2	≥64	≥15
hs-cTnI (Centaur; Siemens)	<3	<8	<7	≥120	≥20
hs-cTnl (Access; Beckman Coulter)	<4	<5	<5	≥50	≥20
hs-cTnl (Clarity; Singulex)	<1	TBD	TBD	≥30	TBD
hs-cTnI (Vitros; Clinical Diagnostics)	<1	TBD	TBD	≥40	TBD
hs-cTnl (Pathfast; LSI Medience)	<3	TBD	TBD	≥90	TBD
hs-cTnl (TriageTrue; Quidel)	<4	TBD	TBD	≥60	TBD



2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation

The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)

Authors/Task Force Members: Jean-Philippe Collet ® * (Chairperson) (France), Holger Thiele ® * (Chairperson) (Germany), Emanuele Barbato (Italy), Olivier Barthélémy (France), Johann Bauersachs (Germany), Deepak L. Bhatt (United States of America), Paul Dendale (Belgium), Maria Dorobantu (Romania), Thor Edvardsen (Norway), Thierry Folliguet (France), Chris P. Gale (United Kingdom), Martine Gilard (France), Alexander Jobs (Germany), Peter Jüni (Canada), Ekaterini Lambrinou (Cyprus), Basil S. Lewis (Israel), Julinda Mehilli (Germany), Emanuele Meliga (Italy), Béla Merkely (Hungary), Christian Mueller (Switzerland), Marco Roffi (Switzerland), Frans H. Rutten (Netherlands), Dirk Sibbing (Germany), George C.M. Siontis (Switzerland)



Population

Patients with suspected MI





W

Vitals 12-lead ECG

Stable patients
No ST-segment elevations


Shock or CPR ST-segment elevations

direct



2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation


The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)





2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation


The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)





2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation

The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)






European Heart Journal (2021) 42, 1289 – 1367 European Society doi:10.1093/eurhearti/ehaa575 **ESC GUIDELINES**

2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation


The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)






2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation

The Task Force for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC)





TROPONIN ASSAYS

TROPONIN ASSAYS



HS troponin assays:

- Superior analytical performance
- Defined by two key criteria:
 - Coefficient of Variation (CV) at the 99th percentile upper reference limit (URL) of <=10%
 - ability to measure concentrations above the limit of detection (LoD) but below the 99th percentile in more than 50% of healthy individuals













HS TROPONIN ASSAYS

-

- Earlier and more accurate diagnosis of AMI¹
- Lower 1-year mortality²
- Efficient risk stratification³
- Operational and Health system advantages
 - o single blood draws, reduced serial testing
 - decrease ED LoS (1.5 3 hour reductions)³
- Patient benefits reduced anxiety, earlier DC
- Personalised medicine capabilities Sex-specific ranges⁴, machine learning tech.
- Cost effective higher test cost, early DC and reduced downstream heath care utilisation⁵












HS Troponin POC assays

What we are currently doing



Hours





POC DEVICES WITH HS ASSAYS





BENCHTOP PLATFORMS

- PATHFAST
- SPINCHIP
- PYLON

PORTABLE HAND HELD PLATFORMS



i-STAT ALINITY

ATELLICA VTLI











Platform	LoD (ng/L)	10%CV (ng/L)	URL Male	URL female	Detectable % healthy reference population	Time to results	Clinical performance studies		
i-STAT Alinity (Abbott)	1.6 (whole blood)	6.9	28	13	11.5	Dr. Cullet	n's STAT		
TriageTrue (Quidel/Ortho)	1.6 (whole blood) 6.9 28 13 1.6 (plasma) 1.9 (whole blood) 8.4 (plasma), 6.2 (whole blood) 6.2 (whole blood) 6.2 (whole blood) 6.2 (whole blood) 6.3 (whole blood) 6.4 (whole blood) 6.5 (whole blood) 6.5 (whole blood) 6.5 (whole blood) 6.6 (whole blood) 6.6 (whole blood) 6.6 (whole blood) 6.7 (whole blood) 6.8 (plasma), 6.9 (plasma), 6.9 (plasma), 6.2 (whole blood) 6.2 (whole blood) 6.2 (whole blood) 6.3 (plasma), 6.4 (plasma), 6.5 (plasma) 6.5 (plasma) 6.5 (plasma) 6.6 (plasma) 6.7 (plasma) 6.8 (plasma), 6.9 (plasma) 6.9 (plasma) 6.9 (plasma) 6.2 (plasma) 6.3 (plasma) 6.2 (plasma) 6.3 (plasma) 6.3 (plasma) 6.4 (plasma) 6.5 (plasma) 6.								
Atellica VTLi (Siemens)	nl is not clear	ed for use on she lia where she	the is a praction of the second secon	cticine i-STAT hi the U.S.	s-Tni is	8min	yes		
i-STAT hs-1	ione in Auseland land land land land land land land	on the on i-S		21	>52%	<17min	yes		
(Spin	ood)	3.7	36.9	27.3	>62%	~10min	yes		
(ET Healthcare)	1.2-1.4	10 (whole blood) 5 (plasma)	27	21	>89%	<20min	no		

	Platform	LoD (ng/L)	10%CV (ng/L)	URL Male	URL female	Detectable % healthy reference population	Time to results	Clinical performance studies
	i-STAT Alinity (Abbott)	1.6 (whole blood)	6.9	28	13	>50%	~15min	yes
	TriageTrue (Quidel/Ortho)	1.6 (plasma) 1.9 (whole blood)	8.4 (plasma), 6.2 (whole blood)	25.7	14.4	>=50%	<20min	yes
	Atellica VTLi (Siemens) 1.2 (plasma) 1.6 (whole blood)		6.7 (plasma) 8.9 (whole blood)	27	18	>=80%	~8min	yes
	Pathfast (LSI Medience)	2.3	15	30	21	>52%	<17min	yes
	SpinChip (SpinChip Diagnostics)	1. (plasma) 2. (whole blood)	3.7	36.9	27.3	>62%	~10min	yes
	Pylon (ET Healthcare)	1.2-1.4	10 (whole blood) 5 (plasma)	27	21	>89%	<20min	no

Study	Assay	Study type	Sample type	
Sorensen NA et al. Clin Chem 2019; 65:1592–601	PATHFAST; LSI Medience	Retrospective derivation	Stored plasma	
Boeddinghaus J, et al. JACC 2020; 75:1111–24; Dakshi AA, et al. EHJ 2023;44.	TriageTrue; Quidel/Ortho Diagnostics	Retrospective derivation Prospective validation	Stored plasma Whole blood (WB)	
Apple FS, et al. Circulation 2022; 146:1918–29; Cullen L, et al. EHJ 2024; 45 (28):2508-2515	Atellica VTLi; Siemens	Prospective derivation/validation	Fresh WB/Stored plasma	
*Greenslade JH, et al. EHJ. ACC. 2025. doi: 10.1093/ehjacc/zuaf068. Epub.	i-STAT; Abbott	Retrospective derivation/validation	Stored plasma	
Koechlin L, et al. JACC 2024;84:726-740.	SpinChip; SpinChip Diagnostics	Retrospective derivation/validation	Stored plasma	

^{*}Greenslade JH, et al relates to off-label use of the i-STAT hs-tnl cartridge. For on-label use please see the i-STAT hs-Tnl cartridge IFU available at https://www.globalpointofcare.abbott/content/dam/ardx/globalpointofcare/apoc/support/i-stat-1/cti-ifu/english-uk/ifu/798049-01A.pdf

JOURNAL OF THE AMERICAN COLLEGE OF CARDIOLOGY

© 2020 THE AUTHORS. PUBLISHED BY ELSEVIER ON BEHALF OF THE AMERICAN

COLLEGE OF CARDIOLOGY FOUNDATION. THIS IS AN OPEN ACCESS ARTICLE UNDER

THE CC BY-NC-ND LICENSE (http://creativecommons.org/licenses/by-nc-nd/4.0/).

Early Diagnosis of Myocardial Infarction With Point-of-Care High-Sensitivity Cardiac Troponin I




Jasper Boeddinghaus, MD, a,b,* Thomas Nestelberger, MD, a,b,* Luca Koechlin, MD, a,b,c Desiree Wussler, MD, a,b,d Pedro Lopez-Ayala, MD, a,b Joan Elias Walter, MD, a,b,d Valentina Troester, MD, a,b Paul David Ratmann, MD, a,b Funda Seidel, MD, a,b Tobias Zimmermann, MD, a,b,d Patrick Badertscher, MD, a,b,e Karin Wildi, MD, a,b,f Maria Rubini Giménez, MD, a,b,g Eliska Potlukova, MD, a,b,d Ivo Strebel, MSc, a,b Michael Freese, RN, a,b Òscar Miró, MD, b,f F. Javier Martin-Sanchez, MD, b,i Damian Kawecki, MD, b,j Dagmar I. Keller, MD,k Danielle M. Gualandro, MD, a,b,d Michael Christ, MD, b,m Raphael Twerenbold, MD, a,b Christian Mueller, MD, a,b for the APACE Investigators

Primary Objective:

• Directly compare the diagnostic accuracy of POC-hs-cTnl-TriageTrue versus two central laboratory assays (Roche Elecsys hs-cTnT and Abbott Architect hs-cTnI)

Secondary Objective:


• The derivation and validation of a POC-hs-cTnl-TriageTrue specific 0/1-h algorithm




The diagnostic accuracy of TriageTrue was high.


Comparable to that provided by the Roche Elecsys hs-cTnT and Abbott Architect hs-cTnI assays

0-hour (presentation) samples



0-hour (presentation) samples






^{*}If chest pain onset >3 h before presentation to the ED

0/1 hour algorithm

Results were comparable to Roche Elecsys hs-cTnT and Abbott Architect hs-cTnI

Safety of early discharge from ED



Circulation

Volume 146, Issue 25, 20 December 2022; Pages 1918-1929 https://doi.org/10.1161/CIRCULATIONAHA.122.061148



ORIGINAL RESEARCH ARTICLE

Single High-Sensitivity Point-of-Care Whole-Blood Cardiac Troponin I Measurement to Rule Out Acute Myocardial Infarction at Low Risk

Fred S. Apple, PhD (i), Stephen W. Smith, MD (i), Jaimi H. Greenslade, PhD, Yader Sandoval, MD (ii), William Parsonage, DM, Isuru Ranasinghe, MBChB, MMed, PhD, Niranjan Gaikwad, MD, MMed, PhD, Karen Schulz, DC, Laura Stephensen, RN, Christian W. Schmidt, MS, Brynn Okeson, MS (ii), and Louise Cullen, MBBS

SAFE EMERGENCY DEPARTMENT DISCHARGE RATE [SEIGE] USA



- 1086 PATIENTS
 - 0 8.1%MI
- WHOLE BLOOD
- <4NG/L PROVIDED:
 - SENSITIVITY OF 98.9% (93.8-100%)
 - NPV OF 99.5% (95% CI: 97.2-100%)
 FOR RULING OUT MI.
 - 30-DAY ADVERSE CARDIAC EVENTS
 0.1% (N=1)





- 1486 PATIENTS
 - o 5.5% MI
- PLASMA



- SENSITIVITY OF 98.8% (93.3-100%)
- NPV OF 99.8% (99.1-100%)
 FOR RULING OUT MI.
 - 30-DAY ADVERSE CARDIAC EVENTS
 0.8% (N=5)







POC hs-cTnl ng/L	<2	3	<4	<5	<6	<7	<8	<9	<10
NPV	100	100	99.8	99.8	99.7	99.5	99.2	99	98.8
	(95.7-100)	(98.8-100)	(99.1-100)	(99.2-100)	(99.1-99.9)	(98.9-99.9)	(98.6-99.6)	(98.3-99.5)	(98-99.3)
Sensitivity	100	100	98.8	97.5	96.3	93.8	88.9	85.2	81.5
	(95.5-100)	(95.5-100)	(93.3-100)	(91.4-99.7)	(89.6-99.2)	(86.2-98)	(80-94.8)	(75.6-92.1)	(71.3-89.2)



CLINICAL RESEARCH

Acute cardiovascular care

Point-of-care high-sensitivity cardiac troponin in suspected acute myocardial infarction assessed at baseline and 2 h

Louise Cullen (1,2,3*, Jaimi Greenslade (1,2,3), William Parsonage (2,4), Laura Stephensen (1,5), Stephen W. Smith⁵, Yader Sandoval⁶, Isuru Ranasinghe^{3,7}, Niranjan Gaikwad⁷, Maryam Khorramshahi Bayat (1,5), Ehsan Mahmoodi⁷, Karen Schulz⁸, Martin Than (1,5), and Fred S. Apple (1,5), on behalf of the SAMIE and SEIGE investigators[†]

Emergency and Trauma Contre, Royal Brisbane and Women's Hoopital, Butterfield Street, Herston, 4019 Queensland, Australia: "Australia Centre for Health Services Innovation, Centre or Health Care from the Individual Centre of Technology, Music Avenue, Kelvin Grove, 4059 Queensland, Australia: "Fiscally of Health Care from Health Care from the Individual Centre, University of Queensland, Herston Road, Herston, 4006 Queensland, Australia: "Department of Cardiology, Royal Brisbane and Women's Hoopital, Brisbane, Australia: "Department of Emergency Medicine at Hennepin Health-care/Hennepin County Medical Center, University of Minnesota School of Medicine, Minneapolis, MN, USA: "Minneapolis Heart Institute, Abbott Northwestern Hospital And Center for Coronary Artery Disease, Minneapolis Heart Institute Foundation, Minneapolis, MN, USA: "Department of Cardiology, The Prince Charles Hospital, Brisbane, Australia: "Department of Laboratory Medicine and Pathology, University of Minnesota, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Christchurch, New Zealand; and ¹⁰Clinical and Forensic Toxicology Laboratory, Hennepin Healthcare/Hennepin County Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Christchurch, New Zealand; and Forensic Toxicology Laboratory, Hennepin Healthcare/Hennepin County Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Christchurch, New Zealand; and Forensic Toxicology Laboratory, Hennepin Healthcare/Hennepin County Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Part Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Part Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Part Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Part Medical Center, Minneapolis, MN, USA: "Emergency Department, Christchurch Hospital, Part Medical Center, Minneapolis, MN, USA: "Emergency Departme

Received 16 October 2023; revised 15 April 2024; accepted 16 May 2024; online publish-ahead-of-print 6 June 2024

See the editorial comment for this article 'Accelerated high sensitivity troponin diagnostics: ready for an even faster pace?', by E. Giannitsis et al., https://doi.org/10.1093/eurheartj/ehae344.

Δ	n	e	TI	ra	•	•

Background and Aims	Strategies to assess patients with suspected acute myocardial infarction (AMI) using a point-of-care (POC) high-sensitivity cardiac troponin I (hs-cTnI) assay may expedite emergency care. A 2-h POC hs-cTnI strategy for emergency patients with suspected AMI was derived and validated.
Methods	In two international, multi-centre, prospective, observational studies of adult emergency patients (1486 derivation cohort and 1796 validation cohort) with suspected AMI, hs-cTnl (Siemens Atellica [®] VTLi) was measured at admission and 2 h later. Adjudicated final diagnoses utilized the hs-cTn assay in clinical use. A risk stratification algorithm was derived and validated. The primary diagnostic outcome was index AMI (Types 1 and 2). The primary safety outcome was 30-day major adverse cardiac events incorporating AMI and cardiac death.
Results	Overall, 81 (5.5%) and 88 (4.9%) patients in the derivation and validation cohorts, respectively, had AMI. The 2-h algorithm defined 66.1% as low risk with a sensitivity of 98.8% [95% confidence interval (CI) 89.3%—99.9%] and a negative predictive value of 99.9 (95% CI 99.2%—100%) for index AMI in the derivation cohort. In the validation cohort, 53.3% were low risk with a sensitivity of 98.9% (95% CI 92.4%—99.8%) and a negative predictive value of 99.9% (95% CI 99.3%—100%) for index AMI. The high-risk metrics identified 5.4% of patients with a specificity of 98.5% (95% CI 96.6%—99.4%) and a positive predictive value of 74.5% (95% CI 62.7%—83.6%) for index AMI.

Primary Objective:

 Derive and validate a 2-h POC hs-cTnl strategy for assessment of patients with suspected AMI.

SAMIE and SEIGE cohorts

POC-hs-cTnl Siemens Atellica VTLi

POC hs-cTnl assessment in suspected AMI 1486 patients with suspected AMI (SAMIE derivation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h ∆ ≥15ng/L Intermediate risk Low risk High risk 5.4% (95% CI 3.4-8.6%) 66.1% (95% CI 61.0-70.8%) 28.6% (95% CI 25.9-31.3%) PPV 74.5% (95% CI 62.7-83.6%) NPV 99.9% (95% CI 99.2-100%) 4.7% (95% CI 3.2-6.8%) with index AMI NPV 98.5% (95% CI 97.7-99.0%) PPV 15.9% (95% CI 10.1-24.2%) 1796 patients with suspected AMI (SEIGE validation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h Δ ≥15ng/L 3 Low risk Intermediate risk High risk 6.1% (95% CI 5.1-7.4%) 53.3% (95% CI 50.9-55.7%) 40.7% (95% CI 38.4-43.0%) NPV 99.9% (95% CI 99.3-100%) 4.5% (95% CI 3.2-6.3%) PPV 49.1% (95% CI 39.8-58.5%) with index AMI NPV 98.0% (95% CI 97.2-98.6%) PPV 10.4% (95% CI 8.5-12.6%)

POC hs-cTnl assessment in suspected AMI 1486 patients with suspected AMI (SAMIE derivation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h ∆ ≥15ng/L Intermediate risk High risk 0 Low risk 5.4% (95% CI 3.4-8.6%) 66.1% (95% CI 61.0-70.8%) 28.6% (95% CI 25.9-31.3%) PPV 74.5% (95% CI 62.7-83.6%) 4.7% (95% CI 3.2-6.8%) INPV 99.9% (95% CI 99.2-100%) with index AMI NPV 98.5% (95% CI 97.7-99.0%) PPV 15.9% (95% CI 10.1-24.2%) 1796 patients with suspected AMI (SEIGE validation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h Δ ≥15ng/L 3 Low risk Intermediate risk High risk 53.3% (95% CI 50.9-55.7%) 6.1% (95% CI 5.1-7.4%) 40.7% (95% CI 38.4-43.0%) NPV 99.9% (95% CI 99.3-100%) PPV 49.1% (95% CI 39.8-58.5%) 4.5% (95% CI 3.2-6.3%) with index AMI NPV 98.0% (95% CI 97.2-98.6%) PPV 10.4% (95% CI 8.5-12.6%)

POC hs-cTnl assessment in suspected AMI 1486 patients with suspected AMI (SAMIE derivation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h ∆ ≥15ng/L Intermediate risk Low risk High risk 5.4% (95% CI 3.4-8.6%) 66.1% (95% CI 61.0-70.8%) 28.6% (95% CI 25.9-31.3%) NPV 99.9% (95% CI 99.2-100%) 4.7% (95% CI 3.2-6.8%) 11 V 77.3% (73% CI 02.7-03.0%) with index AMI NPV 98.5% (95% CI 97.7-99.0%) PPV 15.9% (95% CI 10.1-24.2%) 1796 patients with suspected AMI (SEIGE validation cohort with imputed missing data) 0h <4ng/L 0h ≥60ng/L Other 0h <6ng/L with 2h Δ <5ng/L 2h Δ ≥15ng/L 3 Low risk Intermediate risk High risk 6.1% (95% CI 5.1-7.4%) 53.3% (95% CI 50.9-55.7%) 40.7% (95% CI 38.4-43.0%) NPV 99.9% (95% CI 99.3-100%) PPV 49.1% (95% CI 39.8-58.5%) 4.5% (95% CI 3.2-6.3%) with index AMI NPV 98.0% (95% CI 97.2-98.6%) PPV 10.4% (95% CI 8.5-12.6%)

For consideration

Sample type

Intended users

Real world experience





Improving care by faster risk-stratification by use of next generation point-of-care troponin in patients presenting with possible acute coronary syndrome in the emergency department.

Icare-FASTER trial



doi: 10.1136/bmj.h391 | BMJ 2015;351:h391 | thebmj

BRIEF REPORT

Acute Coronary Syndromes

Emergency department use of a high-sensitivity point-of-care troponin assay reduces length of stay: an implementation study preliminary report

John W. Pickering © ^{1,2}, Laura R. Joyce © ³, Christopher M. Florkowski⁴, Vanessa Buchan⁴, Laura Hamill © ⁵, and Martin P. Than © ^{1,2,6}*

Received 19 June 2024; revised 12 August 2024; accepted 8 October 2024; online publish-ahead-of-print 15 October 2024

¹Department of Medicine, University of Otago Christchurch, Christchurch, New Zealand; ²Emergency Department, Christchurch Hospital, Private Bag 4710, Christchurch 8140, New Zealand; ³Department of Surgery and Critical Care, University of Otago Christchurch, Christchurch, New Zealand; ⁴Canterbury Health Laboratories, Christchurch Hospital, Christchurch, New Zealand; ⁵24hr Surgery, Pegasus Health (Charitable), 401 Madras Street, Christchurch, New Zealand; and ⁶Department of Emergency Medicine, University of Kansas Medical Center, 4000 Cambridge Street Mailstop 1019, Kansas City, KS 66160, USA




Figure 1 Daily mean lengths of ED stay of patients receiving a troponin test for the eight weeks before and after the deployment of the POC test in the ED.

Mean reduction in LoS of **32 mins.**

21.4 fewer patient-hours each day (1196 hours over 8 weeks)

8,000 with suspected ACS

8,000 with suspected ACS

Saving of 32 mins per patient

8,000 with suspected ACS

Saving of 32 mins per patient

Saving of 4,266 hour per year

8,000 with suspected ACS

Saving of 32 mins per patient

Saving of 4,266 hour per year

Saving of 178 days per year



IT

(Codes, LIS, Middleware, results reporting, firewalls)

Equipment

(Purchase orders, distribution, allocation, documentation)

Reagents

(Cartridges, QC material, procurement)

Documentation

(Reports,
Reference change,
work instructions,
trouble shooting,
training records, logs)

Training

(Competency, superusers, forms)

Accreditation


(NATA etc.)


Other


(Bar code scanners, Lock out systems)


• • • •


















Laboratory expertise

Clinical expertise

Project team




Laboratory expertise

Clinical expertise

Project team

Executive leadership








The magic is not in the box.

Practical aspects











"Wait for the formal results."

"That's just POC, we can't act on that!"

"Don't have all the information needed."




Analytical performance

Device performance

Patient needs

ED needs





High-Sensitivity
Cardiac Troponin
Assays:
From Implementation
to Resource Utilization
and Cost Effectiveness

Fig. 1. Practical tips for implementing a hs-cTn assay.

Lee et al. JALM 2024

FAST FORWARD


Opportunities for POC hscTn Assays

FAST FORWARD

Opportunities for POC hscTn Assays

Need evidence for impact of POC hs Troponin assays in many areas of health care delivery





Correct disposition of patients with suspected ACS and + troponin

Primary care


Single RO AMI testing of patients with suspected ACS




Specialist rooms and Outpatients

Assessment and monitoring of patients with acute and chronic conditions

Questions?





Thank you for listening.