

MAYO
CLINIC



Rapid HIV Assays

September 29, 2020

Heart of America Point of Care Network

Brad S. Karon, MD, PhD

Professor of Laboratory Medicine and Pathology

Mayo Clinic

Rochester, MN

Outline

- HIV in US and other countries
- HIV markers and testing algorithm
- Need for rapid HIV in hospitals
- Performance of rapid HIV assays
- Mayo evaluations of rapid antibody and antibody/antigen tests
- Conclusions

HIV in US

- ~ 1,200,000 HIV infected in US as of 2018
- ~ 14% unaware of HIV status (down from close to 25% a decade ago)
- 15,280 deaths in 2018 (all causes)
- New HIV diagnoses 2018:
 - 66% male to male sexual contact
 - 24% heterosexual contact
 - 7% injection drug use
 - 4% male to male sexual contact and injection drug use
 - 42% Black/African American, 27% Hispanic/Latinos, 25% Whites


HIV in US

- Screening for HIV antibodies done with combined antigen.antibody EIA test, takes 3-4 hours technical time, confirmation if positive
 - 8-24 hours TAT standard in most labs
- CDC now recommends HIV screening for all adults and adolescents ages 13-64 in healthcare setting (everyone should get tested at least once), once/year for high risk groups
 - How to get testing, counseling done in order to initiate treatment
 - Rapid HIV tests

HIV outside US

- ~ 37 million living with HIV worldwide in 2018
- ~ 1.7 million new infections
- Around half become infected before age 25
 - 2nd leading cause of death among 20-24 yo
- Sub-Saharan Africa accounts for 2/3 new HIV infections

HIV Testing and Care Continuum



HIV outside US

- Multiple challenges in identifying HIV infected
 - Resources for screening
 - Test and sample stability
 - Resources for confirmation of pos screens
 - Lab personnel for testing
- Rapid HIV
 - Fingertstick or oral fluid sample types
 - Simple methods require minimal training
 - No lab equipment, low cost
 - Long shelf-life and RT storage

Markers of HIV Infection





Figure adapted from Laboratory Testing for the Diagnosis of HIV Infection: Updated Recommendations. CDC 2014.

CDC Testing Algorithm


HIV-1/2 antigen/antibody immunoassay★

Updated:
January, 2018



2018 Quick reference guide: Recommended laboratory HIV testing algorithm for serum or plasma specimens. <https://stacks.cdc.gov/view/cdc/50872>



CDC Testing Algorithm



The FDA-approved single-use rapid HIV-1/HIV-2 Ag/Ab immunoassay can be used as the initial assay in the laboratory HIV testing algorithm for serum or plasma. If any instrumented Ag/Ab test is available, it is preferred due to its superior sensitivity for detecting HIV during acute infection. More data needed on whole blood performance.

What about rapid/POC tests?

- November 2002: FDA approved first rapid HIV diagnostic test kit for use in the U.S.
- OraQuick Rapid HIV-1 Antibody Test
 - Less than a drop of whole blood
 - Results interpreted in 20 minutes
 - 99.6% sensitivity; 100% specificity
 - Room temperature storage
 - 4-6 weeks after infection most enough Ab to test +



POC Testing Characteristics

ADVANTAGES

- Reduce patient loss to follow-up
- Increase access to therapy
- Decrease transmission
- Access more people
- Patient preference

DISADVANTAGES

- Higher rates of misdiagnosis (false pos)
- Longer window period
- Ambiguous test results
- Subjective variability in reading

Evolution of POC Tests

Test Category	HIV Screening Tests	Manufacturer	Run Time	Detects IgG	Detects IgM	Whole Blood	Oral Fluid
Ab test	DPP HIV-1/2 Assay	Chembio	10 min WB 25 min OF	×		×	×
	HIV 1/2 STAT PAK	Alere	15 min	×		×	
	INSTI HIV-1/HIV-2 Ab Test	BioLytical	<2 min	×	×	×	
	OraQuick ADVANCE Rapid HIV-1/2 Ab Test	OraSure Technologies	20 min	×	×	×	×
	Reveal G4 Rapid HIV-1 Ab Test*	MedMira	<2 min	×		×	
	SURE CHECK HIV 1/2 Assay	Chembio	15 min	×		×	
	Uni-Gold Recombigen HIV-1/2	Trinity Biotech	10 min	×	×	×	
Ag/Ab test	Determine HIV-1/2 Ag/Ab Combo	Alere	20 min	×	×	×	

Role of POC HIV Testing at Mayo Clinic (use case for hospital-based rapid HIV)

Labor and Delivery

Rapid testing of women whose HIV status is unknown at labor

Occupational Health

Rapid testing for high- or unknown-risk patients after a needlestick, blood, or body fluid exposure has been reported

HIV Clinic

Anonymous and confidential, rapid testing for patients presenting to the clinic

Other use cases: ED (not done at Mayo), school clinics (one site at Mayo), public health screening



Mother-to-Child HIV Transmission

- Risk of HIV transmission is $<1\%$ when:
 - Recommended antiretroviral/obstetric interventions are used in women who know of infection early in pregnancy
 - Risk $\sim 25\%$ without intervention
- In 2000, of the 6000-7000 HIV-infected women who gave birth in the U.S., 40% had not been diagnosed with HIV before L&D
- 2006 *CDC's revised Recommendations for HIV Testing of Adults, Adolescents, and Pregnant Women in Health-Care Settings*
 - HIV screening included in routine panel of prenatal screening tests for all pregnant women

Rapid Testing for Labor and Delivery



- Starting treatment during L&D or providing it to the newborn within hours after birth can reduce transmission by 50%
- Critical to rapidly obtain HIV test results for women in labor to begin treatment as soon as possible
- CDC recommends rapid HIV testing either prior to the onset of labor or immediately post-partum
 - 9-13% HIV transmission rates achieved when intervention begins intrapartum or neonatal periods
- “Every delivery unit needs to have access to an HIV test that can be done rapidly (i.e., in <1 hour) 24 hours a day”

Rapid Testing for Occupational Exposures

- Occupational transmission of HIV to health care workers is rare
 - As of December 31, 2013, 58 confirmed occupational transmissions of HIV and 150 possible transmissions had been reported in the U.S.
- Post-exposure prophylaxis (PEP) regimens are recommended when occupational exposures to HIV occur
 - Treatment should begin within 2 hours (ideally), but no longer than 72 hours of an exposure
- Rapid testing is critical if HIV status of the source patient is unknown
 - A positive rapid HIV test is preliminarily considered a true positive for the purposes of PEP decision-making

Rapid HIV testing US Public Health

- Overall, sens and spec of waived rapid HIV tests very good in public health setting
 - NYS experience very good except oral fluid testing in one NYS site
 - Some studies found increase false positive with non-lab testing personnel
- Variable effectiveness of screening programs
 - Rapid screening effective in getting HIV results to population that otherwise would not be screened
 - Rates of counseling and confirmation testing vary
 - Rates of entry into medical treatment vary from 47-97%
 - Publicly sponsored programs do better than privately sponsored

Assays/Technologies available

- Immunochromatographic (lateral flow) immunoassay
 - HIV Ag applied to line on nitrocellulose strip
 - Blood diluted in buffer, added to well
 - Lateral flow pulls blood past indicator and over to test line (Ag), then to control line (indicator)



Assays/Technologies available

- Immunochromatographic (lateral flow) immunoassay
- Advantages:
 - Fast (10-20 minutes), few steps, minimal sample processing, varied sample types (whole blood, serum, plasma, oral fluid), many CLIA waived
- Disadvantages:
 - Cannot distinguish HIV-1 and HIV-2, need to dilute sample in buffer, interpretation of lines

Performance of rapid HIV assays

	Sensitivity (95% C.I.)	Specificity (95% C.I.)
OraQuick Advance		
- oral fluid	99.3%(98.4-99.7)	99.8%(99.6-99.9)
- whole blood	99.6%(98.5-99.9)	100%(99.7-100)
- plasma	99.6%(98.9-99.8)	99.9%(99.6-99.9)
Unigold Recombigen	100.0%(99.5- 100)	99.8%(98.3-100)

Performance of rapid HIV assays

- Rapid HIV Ab tests use only one (gp41) or two (gp120 and gp41) antigens as targets
- Despite smaller number Ag targets sensitivity compares well to lab EIA
- One limitation patients treated with high dose antiviral agents (anti gp 41 decreased with therapy)
 - Treated patients unlikely to be tested with rapid HIV test

Performance of rapid HIV assays

- Some studies found poorer detection of HIV Ab early in infection with rapid tests
 - Ab EIA longer window period than Ag/Ab or NAT test (22 vs 17 days after infection)
 - Fewer antigen targets in rapid tests may lead to fewer early pos results compared to lab EIA
 - Some data suggests variability among rapid tests (some may detect IGM Ab better)
 - With non-lab users false positives can be a problem

Performance of rapid HIV assays

- Oral fluid testing
- NYC public screening program 2005-08
- Higher rate of false pos Oraquick results with oral fluid compared to WB
- Many false pos seen in one site, no cause determined
- Still within stated 98% specificity
- CDC now warns that oral fluid testing less sensitive and specific than WB or serum/plasma

Impact of new guidelines for HIV testing

- Start with combined antibody/antigen test
- Only one FDA-approved rapid antigen/antibody combo test available
- Alere Determine™ HIV-1/2 Ag/Ab Combo
- Simultaneous detection of HIV-1 p24 Ag and Ab to HIV-1 and HIV-2
 - Capillary whole blood (waived), serum or plasma
- Rapid differentiation of HIV-1 and HIV-2 antibodies (Multispot no longer manufactured)

Mayo Clinic Evaluation-2005

- Study Design
 - Oraquick, Uni-Gold and Multi-spot evaluated
 - 50 blood bank samples (HIV negative)
 - 20 HIV positive samples (viral load positive)
 - 20 cross-reactive samples (Hep A or B, EBV)
 - 10 EIA positive, WB negative samples
 - Background clarity and line intensity graded
 - 0-3 (0 background best, 3 line intensity best)

Mayo Clinic Evaluation-2005

- HIV negative samples
 - Oraquick: 50/50 negative
 - Uni-Gold: 49/50 negative (1 false positive)
 - Multi-spot: 48/50 negative (2 false positives, undifferentiated)
 - Background clarity:
 - At 10 min, all had 0 background except one sample on Uni-Gold (hemolyzed)
 - At 20 min, Uni-Gold and Multi-spot had higher background than Oraquick

Mayo Clinic Evaluation-2005

- Positive samples
- All 3 methods had 20/20 positive
- All samples resulted in 2-3+ line intensity
- All 3 methods have excellent sensitivity

Mayo Clinic Evaluation-2005

- Crossreactive samples
- 1 sample positive by all 3 methods, record reviewed and re-classified as HIV positive based on history HIV infection
- 19/19 crossreactive samples negative by all 3 methods

Mayo Clinic Evaluation-2005

- EIA Positive, WB negative samples
- 10/10 negative by both Oraquick and Uni-Gold methods
- 6/6 negative by Multi-spot (insufficient volume 4 samples)

Mayo Clinic Evaluation-2005

	OraQuick % (95% CI)*	Uni-Gold % (95% CI)	Multispot % (95% CI)
Specificity	100 (95-100)	98.7 (93-100)	97.3 (91-100)
Sensitivity	100 (84-100)	100 (84-100)	100 (84-100)
Positive Predictive Value	100 (84-100)	95.5 (77-100)	91.3 (72-99)
Negative Predictive Value	100 (95-100)	100 (95-100)	100 (95-100)

Implementation decision (2006-7)

- Implement Oraquick rapid HIV testing using EDTA plasma (ease of use, ease of reading)
- Later changed to EDTA whole blood to standardize between waived and non-waived sites performing testing

Evaluation of Determine™ Ag/Ab Combo Test-2018

- Sample Type: compared 6 waste EDTA whole blood and plasma samples
 - Spiked with positive control material
 - Read in duplicate

	EDTA Whole Blood		EDTA Plasma	
	Reader 1	Reader 2	Reader 1	Reader 2
Clear	3	2	5	6
Difficult	3	4	1	0
Not Clear	0	0	0	0



Evaluation of Determine™ Ag/Ab Combo Test


- Precision: 5 replicates x 5 days using 4 different controls
 - HIV-1 reactive, HIV-2 reactive, p24 reactive, nonreactive
 - Two techs blinded to results read strips
 - 100% concordance with control type
 - 100% concordance between techs

Evaluation of Determine™ Ag/Ab Combo Test

- Accuracy:
 1. Spiked EDTA plasma samples (n = 10) with SeroDetect HIV-1/HIV-2 Ag/Ab Combo Verification Panel
 - Used as reference for testing the Determine™ Ag reaction line
 - 100% concordance (8 Ag reactive/2 Ag non-reactive)
 2. Obtained samples (n = 60) from Hepatitis/HIV Serology Lab with HIV testing performed on Geenius HIV 1/2 Supplemental Assay (Bio-Rad)


Accuracy Results

- Ab Accuracy Compared to Reference Method:
 - 30 positive and 30 negative (reference method)
 - 98% concordance (59/60)
 - Discordant sample was Ab positive by reference method and indeterminate on the Determine™ (invalid control)




Accuracy Results


- Ab Accuracy Compared to Reference Method:
 - 30 positive and 30 negative (reference method)
 - 98% concordance (59/60)
 - Discordant sample was Ab positive by reference method and indeterminate on the Determine™ (invalid control)



**x5 and x10
dilutions
with serum
and saline**



x5 serum



x5 saline

Accuracy Results

- Control line is formed from same components as test lines; consumption of those components by a high titer sample can limit control line intensity—
- **High Dose Hook Effect**

Ambiguity in Lines

- Control lines
 - Any visible pink/red color in control area, regardless of intensity, is considered “valid”
 - If there is no pink/red control line in the control area, even in a pink/red line appears in Ab or Ag area, the result is invalid
 - **Test should be repeated**
 - 8/81 (10%) test have weak/no control lines
 - **May result in increase in repeat testing, potential to not release reactive results**

Ambiguity in Lines

- “Ghost” Lines
 - White abnormalities observed after running strips
 - Predominantly found in a single lot of strips, but present in others
 - Potential to obscure result reading
 - Manufacturing defect (rare)



Non-specific Reactivity

Table 5: Alere Determine™ HIV-1/2 Ag/Ab Combo Reactivity with Specimens from Individuals with Unrelated Medical Conditions and Specimens with Potentially Interfering Substances

Specimen Description	Alere Determine™ HIV-1/2 Ag/Ab Combo (# Reactive/Total Tested)		
	Specificity Testing: Unspiked Samples	Sensitivity Testing: HIV-1 Samples (Weak Reactive)	Sensitivity Testing: p24 Antigen Samples
Human T-cell Lymphotropic Virus (HTLV)	0/10	10/10	10/10
Epstein Barr Virus (EBV)	0/20	10/10	NT
Cytomegalovirus (CMV)	0/20	10/10	10/10
Hepatitis C Virus (HCV)	0/30	10/10	10/10
HBsAg	0/8	NT	NT
Herpes Simplex Virus (HSV)	1/55	20/20	20/20
Syphilis	0/20	10/10	10/10
Toxo IgG	1/55	20/20	20/20
Cancer	2/55	20/20	20/20
Alcoholic Cirrhosis	0/10	10/10	10/10
Flu Vaccine	0/10	10/10	9/9
Anti-HBc	0/10	NT	NT
Multiparous Females	0/10	NT	NT
Drugs	0/10	NT	NT
Hospitalized patients	8/560	55/55	56/56
HAMA	2/54	20/20	20/20
RF	4/150	21/21	21/21
Triglycerides*	3/55	21/21	21/21
Hemoglobin**	0/21	21/21	21/21
Bilirubin**	0/21	21/21	21/21
High Serum Protein**	0/21	21/21	21/21

*Naturally occurring specimens containing more than 500 mg/dL

**Specimens artificially created by adding the potentially interfering substance to normal human serum (Lyophilized hemoglobin: 5 mg/mL; Bilirubin: 0.25 mg/mL; Protein: 0.05 g/mL).

Implementation Decision...

- Agreed to hold off on implementing the Alere Determine™ Ag/Ab Combo Test
 - Challenges and issues observed during the study were concerning
- Continue testing with OraQuick ADVANCE Rapid HIV-1/2 Ab Test
 - Automatic reflex to lab-based Ag/Ab immunoassay
- For rapid testing, means window period longer than lab testing (22 vs 11 days)
- Will need to re-evaluate in future as technologies improve

Conclusions

- HIV still public health threat in US
 - 25% infected globally are still undiagnosed
- Rapid HIV tests have excellent sensitivity for the detection of HIV antibodies
- Rapid HIV testing may be useful for hospital obstetrics, clinic/ED settings, or as part of employee exposure protocol
- Evaluation of devices with intended user group important, potential for false positives greater than with use of automated laboratory tests, as with all POC consider usability and performance data